IonomaX

Tipo III

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Type III

Dentine substitute for bases

IonomaX is a cement containing polycarboxylic acid in its liquid, ensuring chemical adherence to tooth structure. Even though it is an acid, the liquid is biologically compatible and innocuous when placed to dentine and pulp surfaces, due to its high molecular weight. The powder is a complex glass system with high fluorides content capable of being released in the medium where it acts.

Properties: It has specific adhesion and if the powder/liquid ratio is maintained, the product of the reaction is slightly soluble, which enables to obtain a cement stable in the mouth. These features, added to its mechanical properties, make it a filling material (thickness over 0.5 mm) and/or as a dentine substitute wherever dentine must be replaced.

Instructions for Use:

Dosage: Shake the bottle of powder before opening. Dispense two measures of the powder to be used (level tablespoon) and pour three drops of liquid with the dropper bottle, on a smooth glass slab or on non-absorbent paper block. (Weight relation: 2g powder for each gram of liquid)

Mixture: Add all the powder into the liquid in one go or at the most in two times with a Teflon or stainless steel spatula of good quality; mixture time should not exceed 45 seconds. The aspect must be bright and sticky (this ensures adhesion to dental structure) After this time mixture must be discarded.

Application in mouth: Take the material to the clean and dry cavity. Do it without

delay and do not use if material has lost the glossy aspect. This is important to ensure adhesion to the dental surface. The total time available for application in mouth is 1'30" to 2'.

Setting time: Maximum 6 minutes.

Note: Adhesion to dentine can be improved previously treating the surface with a tissue conditioner based on polyacrylate acid.

Presentation:

10 g powder container

10 ml liquid container

Measuring spoon

Storing conditions:

Store in a cool and dry place

EXCLUSIVE USE FOR HEALTH PROFESSIONALS AND INSTITUTIONS.

Important: As the mixture of powder and liquid is very sensitive to environmental conditions such as temperature and humidity, the user must determine if the consistency of the mixture is appropriate. The working and setting times correspond to an room temperature of 23° C. The higher the temperature, the shorter the working and setting times. The opposite happens at lower temperatures.

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LABORATORIOS SL S.A.