IonomaX

Tipo I

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IonomaX

Type I

Glass ionomer cement for fixing stiff restorations and orthodontic bands

IonomaX is a cement containing polycarboxylic acid in its liquid, ensuring chemical adherence to tooth structure. Even though it is an acid, the liquid is biologically compatible and innocuous when placed on dentine and pulp surfaces, due to its high molecular weight. The powder is a complex glass system with high fluorides content capable of being released in the medium where it acts.

Properties: It offers specific adhesion; as long as powder/liquid ratio is maintained, the reaction final product has very low solubility which enables to obtain a stable cement in the buccal medium. These characteristics in addition to its mechanical properties make it a convenient material to place stiff restorations and endodontic bands.

Directions for Use:

Dosage: Shake the bottle of powder before opening. Dispense the amount of powder measures to be used (level spoon) and pour two drops of liquid for each powder measure with the dropper on a flat small slab or a wad of non-absorbing paper. (Weight ratio: 1.8 g powder for each gram of liquid).

Mixture: Add all the powder into the liquid in one go with a Teflon or stainless steel spatula of good quality; mixture time should not exceed 45 seconds. The aspect must be bright and sticky (this ensures adhesion to dental structure) After this time mixture must be discarded. Application in mouth: Apply the material

on the surface of the restoration to adhere. Both surfaces should be clean and dry. Press softly and remove greater excess, and then press hard. Take the rest of excess when the mixture is already hard. The total time available for application in mouth is 1'30" to 2'.

Setting time: Maximum 8 minutes.

Note: Adhesion to dentine can be improved previously treating the surface with a tissue conditioner based on polyacrylate acid.

Presentation:

14 g powder container
10 ml liquid container
Measuring spoon

Storing conditions:

Store in a cool and dry place

EXCLUSIVE USE FOR HEALTH PROFESSIONALS AND INSTITUTIONS.

Important: As the mixture of powder and liquid is very sensitive to environmental conditions such as temperature and humidity, the user must determine if the consistency of the mixture is appropriate. The working and setting times correspond to a room temperature of 23°C. The higher the temperature, the shorter the working and setting times. The opposite happens at lower temperatures.

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LABORATORIOS SL S.A.